


# ZOLL Defibrillators



[ZOLL](#) offers many different defibrillators depending on your customers needs. For AEDs ZOLL offers both semi automatic and automatic units. Most customers prefer the semi automatic. All three offer CPR feed-back. The only difference is the semi automatic user needs to push the shock button. The Rseries ALS defibrillator's most common configuration includes the pacing option. The other options BP, SPO2 and ETCO2 are typically not used in the ASC setting. The ASC's patient monitors already have these features.



[AED3](#) is the ZOLL Flagship AED today. The ZOLL AED 3® features enhanced Real CPR Help® with full-color rescue images and an integrated child mode, making it simple to treat both adult and pediatric victims. [Demo](#)



[G5 Power Heart](#)- With 354 joules of high energy and a synchronized shock this is the most powerful AED today design to give the best shock the first time. The Powerheart® G5 AED with Intellisense™ CPR provides corrective feedback to ensure delivery of high-quality compressions. One button dual-language functionality and comprehensive self-tests give rescuers the confidence to act. [Demo](#)



[AEDPlus](#)—The AED Plus® with clinically proven, integrated, real-time CPR feedback and industry leading long-life consumables, sets the benchmark for dependability, durability, performance and readiness. Bystanders will be able to act quickly and confidently in an SCA emergency. [Demo](#)



[Rseries](#) ALS with Pacing- OneStep Simplicity Delivering therapy is easy with the OneStep™ family of electrodes. Monitor, defibrillate, or pace—without the need for a separate three-lead ECG cable—with OneStep Complete electrodes. To take advantage of Real CPR Help, use OneStep Complete, OneStep CPR, or OneStep Pediatric CPR electrodes. And, with a dual circuit design, OneStep electrodes can remain pre-connected, even during a 30-joule shock test. They're available in anterior/anterior or anterior/posterior configurations.

[Demo](#)

<b><u>ZOLL#</u></b>	<b><u>Description</u></b>
various 210, 228	ZOLL AED PLUS - SEMI-AUTOMATIC
various 210,228 (7=auto)	ZOLL AED PLUS - FULLY AUTOMATIC
8511-001101-01	ZOLL AED 3 Semi-Auto
8511-001102-01	ZOLL AED 3 Fully-Auto
8513-001103-01	ZOLL AED 3 BLS
G5S-80C-S	Powerheart G5 SemiAuto
G5A-80C-S	Powerheart G5 FullyAuto
3502000000-1130012	R Series ALS Package NO EXP
3502000000-1110012	R Series ALS Package W/EXP
3502000000-1130013	R Series BLS Package

Contact your local distributor for more details.