

# SUNTECH

## Oscar 2

24-Hour ABPM

*24-Hour Ambulatory Blood Pressure Monitoring System.*

The Oscar 2 sets the standard in 24-hour ambulatory blood pressure monitoring by utilizing proprietary motion-tolerant technology, optimizing patient comfort, giving clinicians practical interpretive analysis and incorporating the latest monitoring advances. This system comes standard with both the patented Orbit™ ABPM cuff and interpretive software AccuWin™ Pro 4 and NOW features options for central blood pressure measurement and Bluetooth-connected True24 patient diary mobile app.

### *- Optional Features -*



- ◆ **SphygmoCor Inside:** SunTech Medical has partnered with AtCor Medical, the Gold Standard in non-invasive measurement of central BP and arterial stiffness, to bring you the first and only ABP Monitor with SphygmoCor inside. The Oscar 2 ABPM system delivers a patient-focused solution with unmet comfort, data reliability and confidence.
- ◆ **Bluetooth with True24:** The True24 patient diary mobile app\* is designed to pair exclusively with the Oscar 2 ABPM device with optional Bluetooth. This easy-to-use app provides a tutorial on ambulatory blood pressure monitoring and prompts the patient to enter information about their activity, posture or symptoms experienced after a BP measurement is taken. This diary information is extremely valuable, further assisting clinicians with interpreting the ambulatory blood pressure study data and making hypertension treatment decisions.

[Brochure](#)

[Competitive Matrix](#)

[Quick Start Guide](#)

[Instructional Videos](#)



# SUNTECH

---

## **Oscar 2 M250 Monitors**

<b>Part Number</b>	<b>Product Description</b>
<b>SUN-99-0133-00</b>	Oscar 2 ABP Monitor - Standard
<b>SUN-99-0133-01</b>	Oscar 2 ABP Monitor - Bluetooth
<b>SUN-99-0133-02</b>	Oscar 2 ABP Monitor - with SphygmoCor inside
<b>SUN-99-0133-03</b>	Oscar 2 ABP Monitor - with SphygmoCor inside - Bluetooth

## **Oscar 2 M250 Accessories**

<b>SUN-27-0144-A1</b>	AccuWin Pro 4 software (Download)
<b>SUN-91-0130-00</b>	Oscar 2 micro USB cable

## **Oscar 2 M250 Cuffs**

<b>SUN-98-0600-01</b>	OPD- Child cuff (12-19 cm) - Quick connect connector (Box of 5 cuffs)
<b>SUN-98-0600-41</b>	OPD- Child cuff (12-19 cm) - Bayonet connector (Box of 5 cuffs)
<b>SUN-98-0239-01</b>	#1 Orbit comfort cuff (18-26 cm) - Quick connect connector
<b>SUN-98-0239-02</b>	#2 Orbit comfort cuff (26-34 cm) - Quick connect connector
<b>SUN-98-0239-03</b>	#3 Orbit comfort cuff (32-44 cm) - Quick connect connector
<b>SUN-98-0239-04</b>	#4 Orbit comfort cuff (42-55 cm) - Quick connect connector
<b>SUN-98-0239-21</b>	#1 Orbit comfort cuff (18-26 cm) - Bayonet connector
<b>SUN-98-0239-22</b>	#2 Orbit comfort cuff (26-34 cm) - Bayonet connector
<b>SUN-98-0239-23</b>	#3 Orbit comfort cuff (32-44 cm) - Bayonet connector
<b>SUN-98-0239-24</b>	#4 Orbit comfort cuff (42-55 cm) - Bayonet connector
<b>SUN-91-0028-22</b>	Patient hose adapter, female bayonet to female quick connect

**For more details contact your local distributor**