

SUNTECH

Tango M2

Stress BP

Automated Blood Pressure Monitor for Cardiac Stress and Exercise Testing

SunTech Medical's Tango M2 automated NIBP monitor uses DKA technology specifically engineered for the cardiac stress testing environment. With advanced proprietary algorithms SunTech Medical stands alone in the market of motion tolerant BP measurement. The Tango M2 allows for the focus to be on the patient by automatically taking accurate and reliable BP measurements rather than spending valuable time trying to take measurements manually. Tango M2 can add automated BP and SpO2 to your stress ECG system, creating a complete cardiac diagnostic testing center.



[Brochure](#)

[User Manual](#)

[Quick Start Guide](#)

[Instructional Videos](#)

[Interface Instructions](#)

- Features and Benefits -

- ◆ **Seamless Integration** Automated communication with your stress system reduces the risk of transcription errors.
- ◆ **Non Exercise BP Mode** Allows BP measurement during patient set-up and recovery without an ECG signal.
- ◆ **Stat Mode** Rapidly repeated automatic BP measurements for time sensitive and emergent situations.
- ◆ **Color LCD** Improved usability with a new 7" color LCD display.
- ◆ **Verify Measurements** See & hear the Korotkoff sounds using the onscreen display and the included headphones.
- ◆ **Data Retrieval** Easier troubleshooting with 300 BP reading history and USB capabilities for measurement retrieval
- ◆ **Field Upgrades** USB port allows for field upgrades, making sure end-users always have the current software.
- ◆ **E-Library** Onboard E-Library, including step-by-step interface notes and tutorials.

- Accessory Options -

Single Use Cuff Kits
Reusable Orbit K Cuffs
Patient Cables
Orbit K Cuff Microphones
Preventative Maintenance Kits



SUNTECH

Tango M2 Base Monitors

Item Number	Product Description
SUN-99-0088-00	Tango M2 Monitor Without ECG
SUN-99-0088-01	Tango M2 Monitor with ECG

Tango M2 System Monitors

SUN-99-0132-00	System for GE case, case 8000
SUN-99-0132-01	System for Cabridge Heart 2000, Cabridge HearTWave II
SUN-99-0132-04	System for Quinton Q-Stress (up to version v4.6), Q-4500
SUN-99-0132-05	System for Mortara, X-Scribe, Philips StressVue, Quinton Q-Stress (version 6+)
SUN-99-0132-06	System for WA Cardioperfect, Midmark Iqmark EZStress
SUN-99-0132-07	System for Philips ST80i

Orbit K- Reusable BP Cuffs

SUN-98-0061-01	Small adult Orbit-K cuff (18-27cm)
SUN-98-0061-02	Adult Orbit-K cuff (25-35 cm)
SUN-98-0061-03	Large adult plus Orbit-K cuff (32-44 cm)
SUN-98-0061-05	Adult plus Orbit-K cuff (27-40 cm)
SUN-98-0062-21	Small adult Orbit-K cuff with microphone (18-27cm)
SUN-98-0062-22	Adult Orbit-K cuff with microphone (25-35 cm)
SUN-98-0062-23	Large adult plus Orbit-K cuff with microphone (32-44 cm)
SUN-98-0062-25	Adult plus Orbit-K cuff with microphone (27-40 cm)

Single Patient Use Cuff Kits (Box of 20 cuffs and mic pads)

SUN-98-0700-01	Small adult (17-25 cm)
SUN-98-0700-02	Adult (23-33 cm)
SUN-98-0700-03	Adult long (23-33 cm)
SUN-98-0700-04	Large adult (31-40 cm)
SUN-98-0700-05	Large adult long (31-40 cm)

Tango M2 Accessories

SUN-99-0027-39	Preventative maintenance kit - (15" patient cable, adult plus, and large adult cuffs)
SUN-99-0127-01	15" patient cable
SUN-99-0235-01	Tango M2 microphone (12")
SUN-99-0235-00	Tango M2 microphone (18")

For more details contact your local distributor