

VIOSCAN®

PORTABLE BLADDER SCANNER



[Demo Video](#)

The Infinium VioScan® is a portable, non-invasive ultrasound bladder scanning device. The VioScan® provides a patented algorithm to measure and determine bladder volume and the need for catheterization. The VioScan® features real time imaging and accurate volume measurement with an easy positioned hand-held probe. A high definition color screen and built in strip chart recorder are standard.

The VioScan® offers cost saving benefits in Urology, Gynecology, General Surgery, ICU, Emergency Medicine, and Nursing Home Facilities.

Infinium **VIOSCAN®** Portable Bladder Scanner:

- 8.4 inch hi-definition color screen
- Data storage of 100 cases
- Rechargeable battery (2 included)
- Optional rolling cart



Item ID	Product Description
INF-VIS906	VIOSCAN PORTABLE BLADDER SCANNER. ONE YEAR WARRANTY
INF-VIS906MS	VIOSCAN BLADDER SCANNER MOBILE STAND
INF-VIS906RP	VIOSCAN BALDDER SCANNER RECORDING PAPER. 10 ROLLS PER BOX

For more information contact your local distributor.