



TelScope® Telehealth System

Intraoral examination and telehealth image and video device
Compatible with all smart devices with TelScope® app



TelScope Telehealth System is the tool and app combination that allows single-handed, illuminated intraoral examinations and turns any smart device into an intraoral camera to capture superior quality intraoral telehealth images and video.


PRODUCT FEATURES

- Places light source inside the mouth via the illuminating, clear depressor, providing the **best view inside the mouth**
- **Frees up an extra hand** during intraoral examination
- **Antimicrobial** handle
- **Rechargeable**, universal handle attaches to any smart device
- **Portable** and lightweight



APP FEATURES

- Programmed to capture medical-grade, **superior quality images and videos** inside the mouth and the throat
- Measure, circle, & comment on areas of concern in the mouth
- **Securely send** encrypted, password-protected files to a dentist or doctor **instantly**
- Available for Android and iOS



Technical Specifications

Dimensions (Handle)	47 x 110 x 165 mm
Dimensions (Depressor)	73 x 200 x 180 mm
Weight (Handle)	0.20 kg
Weight (Depressor)	0.50 kg
File storage (App)	Store to phone or in app; metadata wiped
File transmission (App)	Send securely in password-protected folders
Software	App available for iOS and Android
Warranty	2 years