

ULTRAFLO[®]-R PUSH BUTTON BLOOD COLLECTION SET



1

SAFETY WITH THE PUSH OF A BUTTON

- Single-handed activation retracts the needle instantly from the vein into the needle hub
- Minimizes sharp exposure

2

INCREASED FLOW RATE

- Ultraflo[®]-R (Retract) Needle Cannula increases flow rate 200% (23G) and 60% (21G)*
- Improve patient comfort: Use 23G instead of 21G without compromising flow rate
- Approved for short-term infusion (<2 hrs)

3

BEVEL TECHNOLOGY

- Unique bevel design requires lower penetration force*
- More comfortable patient experience

*Compared to leading-brand standard cannula



INNOVATIVE MED TECH

w: InnovativeMedTech.com | t: 866-560-1820



SAFER & EASIER

Designed to bring comfort for patients and clinicians

REF #	ULTRAFLO [®] -R PUSH BUTTON BLOOD COLLECTION AND INFUSION SET	BOX	CASE QTY
365	21g x 3/4" UF needle with 12" tubing and luer adapter	50	200
364	23g x 3/4" UF needle with 12" tubing and luer adapter	50	200
363	25g x 3/4" needle with 12" tubing and luer adapter	50	200
338	21g x 3/4" UF needle with 7" tubing and luer adapter	50	200
336	23g x 3/4" UF needle with 7" tubing and luer adapter	50	200
335	25g x 3/4" needle with 7" tubing and luer adapter	50	200
326	21g x 3/4" UF needle with 12" tubing w/o luer adapter	50	200
324	23g x 3/4" UF needle with 12" tubing w/o luer adapter	50	200
323	25g x 3/4" needle with 12" tubing w/o luer adapter	50	200
687	21g x 3/4" UF needle with 12" tubing and luer adapter w holder attached	25	200
688	23g x 3/4" UF needle with 12" tubing and luer adapter w holder attached	25	200
689	25g x 3/4" needle with 12" tubing and luer adapter w holder attached	25	200



INNOVATIVE MED TECH

w: InnovativeMedTech.com | t: 866-560-1820